

# De novo spatial proteomic profiling of immune synapses using machine learning-guided Microscoop You-Pi Liu<sup>\*</sup>, Weng Man Chong<sup>\*</sup>, Harry Huang<sup>\*</sup>, Yu-Chih Lin, Yi-De Chen, Chantal Hoi Yin Cheung, Chia-Wen Chung, Chih-Wei Chang and Jung-Chi Liao<sup>#</sup>

### For more information, please visit https://syncell.com/

### Abstract

The spatial proteome of the immune synapse (IS) between a target cell and a lymphocyte is fundamentally important to understand the mechanism of cell-mediated immunity for both immuno-oncology and therapeutic applications. In this research, we used Microscoop<sup>TM</sup>, a fully-automatic microscope system integrated with a machine learning-based algorithm, to best determine the proteome of the IS. We used Raji B cells as antigen-presenting cells (APCs) and induced the formation of IS by incubating with Jurkat T cells. Multiple IS images were applied to train our algorithm using convolution neural network-based deep learning. A sequential process including fluorescence imaging, deep learning-enabled pattern generation, and photochemical labeling was implemented to achieve spatial-specific biotinylation of the IS proteins. Moreover, Microscoop is capable of repeating the process automatically on thousands of fields of view, such that sufficient amount of immune-synaptic proteins are biotinylated for protein identification with mass spectrometry. We have successfully labeled and isolated proteins from spreading T-cell and spatially reorganized interfaces between T cells and APCs. Following MS-based proteomic analysis, several hundreds of proteins were identified, including the proteins known to be specifically associated with T cell receptor (TCR) activation such as LcK, one of the major factors involved in TCR signaling at IS. More interestingly, we identified several proteins novel for IS, including proteins involved in phosphatidylinositol signaling. Our data showcases the capability of subcellular de novo spatial proteomics of the Microscoop technology, revealing novel factors responsible for initiating the immune response of a lymphocyte and shedding light on immune checkpoint signaling and tumor immunotherapy.

## Spreading T cells and T-B cell conjugation models are used to reveal novel biomarkers on IS



## Spatial proteomic profiling via machine learning-guided microscopic photolabeling



SYNCELL Inc., Taipei, Taiwan

Schemes show the cellular localization of IS in a T cell spreading assay (A<sub>i</sub>) and T-B cell conjugation (A<sub>ii</sub>). IS is also known as supramolecular activation cluster, SMAC (A<sub>iii</sub>), This structure is composed of segregated rings with cluster of proteins, including the central ring, cSMAC; the peripheral ring, pSMAC and the distal ring, dSMAC. (B) In the spreading assay, cells were seeded on the CD3 and CD28 antibodies-coated coverslips for 15 min, and fixed and stained with indicated antibodies and dyes. Side view of the cell in the inset in  $B_i$ . IS indicated by CD3 is present at the bottom of the cell. (C) In T-B conjugation, T cells were first seeded on poly-D lysine coated coverslips. After incubating for a few minutes, SEE-activated B cells were added and incubated for another 10 minutes. IS denoted by CD3 staining is indicated by arrows. Bars=  $5 \mu$ M.

raw image

AI defined

merge



Cells or tissue samples were imaged and recognized by computer vision or AI (Step 1); Samples were incubated with photoactivatable probes, and the photolabeling was driven by two-photon illumination (Step 2); Proteins in these labeling regions were tagged with probes and further purified by an affinity pulldown assay (Step 3); De novo biomarkers were identified by LC-MS/MS (Step 4).



Modified from Courtney et al., 2018. doi.org/10.1016/j.tibs.2017.11.008

<sup>#</sup>corresponding author (Email: jcliao@syncell.com) \* These authors contribute equally to this work

## Immune synapses are recognized by AI and precisely labeled by photoreactive probes



Spatial components of immune synapse are identified using Microscoop

labeling on IS ( $E_{i-ii}$ ). Bars= 5  $\mu$ M.



nown
CD3D
CD3E
LcK
ΙΤΚ
Zap70
РІЗК
Csk
NFKB

### De novo LCP2 INPP4A GRAP2 PTPN6 Galectin-9 VAV family

Molecules involved in T cell receptor (TCR) signaling, and a list of potential de novo components are identified using Microscoop.



For more information, please visit our website