



Microscop, a Discovery-Based Image-Guided Proteomics Technology, Reveals Novel Factors on Amyloid- β Aggregates in Differentiated SH-SY5Y Cells

Chien-Chang Huang, Weng Man Chong, Chun-Kai Huang, Yi-De Chen, Chia-Wen Chung, Hsiao-Jen Chang, Jung-Chi Liao*


SYNCELL Inc., Taipei, Taiwan

For more information, please visit <https://syncell.com/> *corresponding author (Email: jliao@syncell.com)


Abstract :

Aggregation of amyloid- β peptides ($A\beta$) is a hallmark of Alzheimer’s disease (AD). Subcellular distribution of $A\beta$ is well recognized under a microscope, but its pathological and physiological functions remain unclear, partially due to our limited understanding of interacting proteins and corresponding signaling pathways associated with $A\beta$. Existing spatial proteomics technologies focus on mapping of known proteins using antibody panels/arrays, hindering hypothesis-free proteome discovery. In this study, we used Microscop™, a fully-automated microscopy-guided subcellular photolabeling with a machine learning-based precision recognition to enable discovery-based image-guided proteomics. We applied $A\beta$ 1-42 deposition in human neuroblastoma, SH-SY5Y, differentiated cells as an AD model. Multiple images of $A\beta$ 1-42 aggregates were applied to segment to locations of specific $A\beta$ 1-42 aggregates of interest using convolutional neural networks (CNN)-based deep learning. Microscop was used to illuminate these segmented regions to induce photochemical reactions of proprietary photosensitive probes and trigger spatial covalent labeling of proteins adjacent to $A\beta$ 1-42 aggregates. This spatial-specific photochemical labeling process was repeated automatically on thousands of microscopic fields of view to accumulate enough $A\beta$ 1-42-associated proteins for LC/MS-MS-based proteome identification. A series of novel factors were discovered to be associated with $A\beta$ aggregation in SH-SY5Y differentiated cells. We further validated these newly identified proteins using antibody staining and found that these proteins indeed colocalized with $A\beta$ 1-42. The finding of these novel factors opens a door to reveal associated signaling pathways related to AD. They may also serve as new diagnostics biomarkers or new AD drug target. Our study not only reveals the $A\beta$ -associated spatial proteome, but also demonstrates its possible broad applications on discovery-based spatial proteomics in neuroscience.

Microscop™ : A novel platform for spatial proteomic profiling via machine learning-guided microscopic photolabeling




Proteomic discovery assay for in vitro amyloid- β plaques




Locations of $A\beta$ 1-42 aggregates were segmented in real time using Microscop

A modified otsu thresholding method was first applied to images after image pre-processing. A following white spot extraction method by logic operation was implemented to segment $A\beta$ 1-42 aggregates with high precision (Fig. A-C).




Targeted photo-induced biotinylation using Microscop

SH-SY5Y cells were differentiated for 5 Days. $A\beta$ 1-42 were expressed by lentivirus infection for 72 hrs. The $A\beta$ 1-42 aggregated spots were indicated by immunostaining. The cells were then incubated with a photo-activatable probe, and the spatial photolabeling process was performed at the segmented locations for thousands of fields of view fully automated with Microscop to biotinylate proteins in the regions of $A\beta$ 1-42 aggregates. Confocal imaging validated the precision of photolabeling (Fig. D-F). Fig. G-I indicates precise labeling viewed from the z direction. High-speed photolabeling for thousands of fields of view was necessary to assure collection of a large number of proteins of interest enough for mass spectrometry sensitivity.




Precise photolabeling at Z-axis



Spatial protein composition of $A\beta$ 1-42 aggregates were identified by Microscop

LC-MS/MS results of well-known $A\beta$ -related proteins



The proteome of the photo-labeled sample by two-photon (2P) illumination was compared with that of the unlabeled sample (CTL). In total, 837 proteins were identified in the 2P enriched group. 24 proteins were found to map to the existing $A\beta$ -related protein database in UniProt. Among them, 15 proteins are well-studied and highly related to $A\beta$ aggregates. Further validation of high-ranked candidates by immunofluorescence staining will be performed to identify possible novel $A\beta$ -associated proteins.

Drebrin	ADNP	VDAC1	3-hydroxyacyl-CoA dehydrogenase type-2
CDK5	Reticulon-3	Tmed10	Delta(24)-sterol reductase
TPK YES	Reticulon-4	Gelsolin	Ubiquitin-60S ribosomal protein L40
TPK LYN	hnRNP A2B1	Histone H2B type 2-E	Amyloid- β precursor protein
Rab11A	PARP1	Histone H2A type 2-C	Tubulin alpha-1B chain
Pin1	OSTC	Histone H4	Clathrin heavy chain 1
SynAz1*			
SynAz2*			

Conclusion

SYNCELL Microscop is a new technology platform for hypothesis-free spatial protein biomarker discovery. It has been used to identify novel stress granule proteins in bulk in another study. Here we have precisely biotinylated proteins in $A\beta$ 1-42 aggregates of the SH-SY5Y differentiated cell model in high content using Microscop. We have further performed LC-MS/MS and identified the proteome of these $A\beta$ 1-42 aggregates. This proteome not only contained known $A\beta$ -associated proteins, but also repeatedly revealed proteins that are not considered as $A\beta$ -related proteins so far. Further validation is in progress to check whether novel $A\beta$ -associated biomarkers can be identified by SYNCELL Microscop.

SYNCELL

